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Referências

- Vignaux L, Grazioli S, Piquilloud L, Bochaton N, Karam O, Jaecklin T, Levy-jamet Y, Tourneux P, Jolliet P, Rimensberger P. Optimizing patient ventilator synchrony during invasive ventilator assist in children and infants remains a difficult task. *Ped Crit Care Med*. 2013;14(7), 316-325.
- Beck J, Reilly M, Grasselli G, et al. Patient-ventilator interaction during neurally adjusted ventilatory assist in low birth weight infants. *Pediatr Res* 2009;65(6):663-8.
- Longhini F, Ferrero F, De Luca D, et al. Neurally adjusted ventilatory assist in preterm neonates with acute respiratory failure. *Neonatal*. 2015;107(1):60-7
- Mally PV, Beck J, Sinderby C, et al. Neural breathing pattern and patient-ventilator interaction during neurally adjusted ventilatory assist and conventional ventilation in newborns. *Pediatr Crit Care Med* 2018;19(1):48-55.
- Beck J, Emeriaud G, Liu Y, Sinderby C. Neurally-adjusted ventilatory assist (NAVA) in children: a systematic review. *Minerva Anestesiol* 2016;82(8):874-83.
- Beck J, Sinderby C. Neurally adjusted ventilatory assist in newborns. *Clin Perinatol* 2021, Dez;48(4):783-811.
- Tabacaru CR, Moores Jr RR, Khoury J, Rozycki HJ. NAVA-synchronized compared to nonsynchronized noninvasive ventilation for apnea, bradycardia, and desaturation events in VLBW infants. *Pediatr Pulmonol*. 2019 Nov;54(11):1742-6
- Hovespyan K, Firestone KS, Moore J, Stein H. Effects of NAVA Compared to SIMV Ventilation on Cardiac Function in Preterm Neonates. *Resp Care* 2020;65(10):3451-491.
- Surkov D M. Neurally Adjusted Ventilatory Assist Mode of Mechanical Ventilation in Neonates with Hypoxic-Ischemic Encephalopathy. *Ukr J Cardio Surg* 2019;2(35):102-11.
- Baudin F, Emeriaud G, Essouri S, Beck J, Javouhey E, Guerin C. Neurally adjusted ventilatory assist decreases work of breathing during non-invasive ventilation in infants with severe bronchiolitis. *Crit Care* 2019;23(1):120.
- De la Oliva P, Schuffelmann C, Gomez-Zamora A, Vilar J, Kacmarek RM. Asynchrony, neural drive, ventilatory variability and COMFORT: NAVA vs pressure support in pediatric patients. A randomized cross-over trial. *Int Care Med*. 2012;38:838-46.
- Piastra M, De Luca D, Costa R, et al. Neurally adjusted ventilatory assist vs pressure support ventilation in infants recovering from severe acute respiratory distress syndrome: nested study. *J Crit Care* 2014;29(2):312:e1-5.
- Firestone KS, Beck J, Stein H. Neurally Adjusted Ventilatory Assist for Non-invasive Support in Neonates. *Clin Perinatol*. 2016 Dez;43(4):707-24.
- Chidini G, De Luca D, Calderini E, et al. Implementation of noninvasive neurally adjusted ventilatory assist in pediatric acute respiratory failure: a controlled before-after quality improvement study. *J Anesth Analg & Crit Care* 2021(Set);1:1-9
- Kallio M, Peltoniemi O, Anttila E, Pokka T, Kontiokari T. Neurally Adjusted Ventilatory Assist (NAVA) in Pediatric Intensive Care – A Randomized Controlled Trial. *Pediatr Pulmonol*. Jan/2015; 50(1):55-62.
- Lee J, Kim HS, Jung YH, et al. Neurally adjusted ventilatory assist for infants under prolonged ventilation. *Pediatr Int* 2017;59(5):540-4
- Baez Hernandez N, Milad A, Li Y, et al. Utilization of neurally adjusted ventilatory assist (NAVA) mode in infants and children undergoing congenital heart surgery: a retrospective review. *Pediatr Cardiol* 2019;40(3):563-9.
- Gibu CK, Cheng PY, Ward RJ, et al. Feasibility and physiological effects of noninvasive neurally adjusted ventilatory assist in preterm infants. *Pediatr Res* 2017; 82(4):650-7.
- Parikka V, Beck J, Zhai Q, Leppäsalo J, Lehtonen L, Soukka H. The effect of caffeine citrate on neural breathing pattern in preterm infants. *Early Hum Dev*. Out/2015;91(10):565-8.
- Rahmani A, Ur Rehman N, Chedid F. Neurally adjusted ventilatory assist (NAVA) mode as an adjunct diagnostic tool in congenital central hypoventilation syndrome. *J Coll Physicians Surg Pak* 2013; Fev;23(2):154-156
- Makker K, Cortez J, Jha K, et al. Comparison of extubation success using non-invasive positive pressure ventilation (NIPPV) versus noninvasive neurally adjusted ventilatory assist (NI-NAVA). *J Perinatol* 2020;40(8):1202-10.
- Keszler M. Volume-targeted ventilation: one size does not fit all. Evidence based recommendations for successful use. *Arch Dis Child Fetal Neonatal Ed* 2018;1-5
- Lee BK, Shin SH, Jung YH, et al. Comparison of NIV-NAVA and NCPAP in facilitating extubation for very preterm infants. *BMC Pediatr* 2019 Ago 28;19(1):298
- Morita et al. The usability of ventilators: a comparative evaluation of use safety and user experience *Crit Care* (2016) 20:263
- Ducharme-Crevier L, Du Pont-Thibodeau G, Emeriaud G. Interest of Monitoring Diaphragmatic Electrical Activity in the Pediatric Intensive Care Unit. *Crit Care Res Pract*. 2013: 384210.
- Stein H, Firestone K. Application of neurally adjusted ventilatory assist in neonates. *Semin Fetal Neonatal Med*. Fev/2014;19(1):60-9.
- Emeriaud G, Larouche A, Ducharme-Crevier L et al. Evolution of inspiratory diaphragm activity in children over the course of the PICU stay. *Intensive Care Med*. 2014 Nov;40(11):1718-26
- Protain AP, Firestone KS, McNinch NL, Stein HM. Evaluating peak inspiratory pressures and tidal volume in premature neonates on NAVA ventilation. *Eur J Pediatr*. Jan/2021;180(1):167-175
- Parikka V, Beck J, Zhai Q, Leppäsalo J, Lehtonen L, Soukka H. The effect of caffeine citrate on neural breathing pattern in preterm infants. *Early Hum Dev*. Out/2015;91(10):565-8
- Amigoni A, Rizzi G, Divisic A et al. Effects of propofol on diaphragmatic electrical activity in mechanically ventilated pediatric patients. *Intensive Care Med*. Out/2015;41(10):1860-1.
- Lee J, Parikka V, Lehtonen L, Soukka H. Parent-infant skin-to-skin contact reduces the electrical activity of the diaphragm and stabilizes respiratory function in preterm infants. *Junho/2021 Pediatr Res*.
- Baudin F, Emeriaud G, Essouri S et al. Physiological effect of prone position in children with severe bronchiolitis: a randomized cross-over study (BRON-CHIO-DV). *J Pediatr* 2019;205:112-9.
- Wolf G, Walsh B, Green M, Arnold J. Electrical activity of the diaphragm during extubation readiness testing in critically ill children. *Pediatr Crit Care Med* 2010;12:e220e4.
- Poets CF, Roberts RS, Schmidt B, et al. Canadian Oxygen Trial Investigators. Association between Intermittent Hypoxemia or Bradycardia and Late Death or Disability in Extremely Preterm Infants. *JAMA*. 11 Ago/2015;314(6):595-603
- Morley CJ, Davis PG, Doyle LW et al. COIN Trial Investigators. Nasal CPAP or intubation at birth for very preterm infants. *N Engl J Med*. 14 Fev/2008;358(7):700-8.
- SUPPORT Study Group of the Eunice Kennedy Shriver NICHD Neonatal Research Network, Finer NN, et al. Early CPAP versus surfactant in extremely preterm infants. *N Engl J Med*. 27 Mai/2010;362(21):1970-9.
- Courtney SE, Pyon KH, Saslow JG, et al. Lung recruitment and breathing pattern during variable versus continuous flow nasal continuous positive airway pressure in premature infants: an evaluation of three devices. *Pediatrics*. Fev/2001;107(2):304-8.
- Dargaville PA, Gerber A, Johansson S, et al. Incidence and Outcome of CPAP Failure in Preterm Infants. *Pediatrics*. 2016;138:e20153985-e20153985.
- Houtekie L, Moerman D, Bourleau A, et al. Feasibility Study on Neurally Adjusted Ventilatory Assist in Noninvasive Ventilation After Cardiac Surgery in Infants. *Respir Care*. Jul/2015;60(7):1007-14
- Longhini F, Scarlino S, Gallina MR et al. Comparison of neurally-adjusted ventilator assist in infants before and after extubation. *Minerva Pediatr* Abr/2018;70(2):133-40
- Yagui ACZ, Goncalves PA, Murakami SH, et al. Is noninvasive neurally adjusted ventilatory assistance (NIV-NAVA) an alternative to NCPAP in preventing extubation failure in preterm infants? *J Matern Fetal Neonatal Med* 2019;1-151
- Shetty S, Evans, K, Cornuauud P, et al. Neurally Adjusted Ventilatory Assist in very prematurely born infants with evolving/established bronchopulmonary dysplasia. *AJP Rep* 2021;11(4):e127-e131.
- Lehtonen, L. (EPNV, 2014). NAVA experiences and research in preterm infants. Extraído de <http://www.criticalcarenews.com>.
- Lehtonen, L. (EPNV, 2014). Hospital in Finland experiences a weight gain of 35% with NAVA - neonatal NAVA and individualizing treatment at bedside. Extraído de <http://www.criticalcarenews.com>
- Cosi G, Monzani A, Genoni G et al. Neurally Adjusted Ventilatory Assist (NAVA): a prospective interventional study in neonates. *Minerva Pediatr*. 05 abril 2019
- Stein H, Howard D. Neurally Adjusted Ventilatory Assist (NAVA) in Neonates less than 1500 grams: a retrospective analysis. *J Pediatr* 2012;160:786e9.
- Meyers, N, Rodrigues N, Ari A. High-frequency oscillatory ventilation: A narra-tive review. *Can J Resp Ther* 2019. 55:40-6
- Bordessoule A, Piquilloud L, Lyazidi A, Moreira A, Rimensberger PC. Imposed Work of Breathing During High-Frequency Oscillatory Ventilation in Spontaneously Breathing Neonatal and Pediatric Models. *Resp Care Set/2018*, 63(9):1085-1093.
- Stein H. (APA, 2014). Neonatal outcomes. Extraído de www.criticalcare-news.com.
- Galindo-Filho, V.C. et al. Radioaerosol Pulmonary Deposition Using Mesh and Jet Nebulizers During Noninvasive Ventilation in Healthy Subjects. *Respir. Care* 2015, 60(9):1238-1246

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Fabricante · Maquet Critical Care AB · Röntgenvägen 2 SE-171 54 Solna · Suécia · +46 (0)10 335 73 00